



Chain Drive Sash Option

Air Master Systems introduces a superior Chain Driven Sash option for Millennia Fume Hoods

- Adaptable to all MAF/MFF/MFS/MFD Fume hoods.
- Will be standard on 8' (96") wide Millennia bench hoods.
- Provides a smooth, easy, one handed sash operation from extreme right hand or left hand side.
- Value Priced.
- Uses a full length bar, front and rear, for stability and longevity.
- Uses metal sprockets and metal chain for durability.
- ADA compliant.
- Delivery from receipt of order is our standard 4 weeks.
- Meets all SEFA standards and codes.

© AIR MASTER SYSTEMS CORP.
6480 Norton Center Dr. • Muskegon, MI 49441
Phone 231.798.1111 • Fax 231.798.4000

On the Web at WWW.AIRMASTERSYSTEMS.COM



Chain Drive Sash Option



The Chain and Sprocket driven sash can be raised and lowered with minimal effort with one hand, from the right or left side of the sash.

For the months of August and September 2006 all purchase orders received for this product will be priced at our base pulley/cable price.

Direct all inquiries to: Air Master Systems Corporation,
6480 Norton Center Drive,
Muskegon, MI 49441
Phone 231-798-1111 or Fax 231-798-4000
Attention - Estimating or Customer Service

© AIR MASTER SYSTEMS CORP.
6480 Norton Center Dr. • Muskegon, MI 49441
Phone 231.798.1111 • Fax 231.798.4000

On the Web at WWW.AIRMASTERSYSTEMS.COM